

VET – 500mA Generator – 400FS Techniques – Table 2 Pg 1

Region Name: Pelvis									
FEL/VD		SID 40		Table Grid		Lg Focal		AEC: Off	
cm	6	8	10	12	14	16	18	20	22
kVp	50	52	56	60	63	67	71	75	80
mAs	7.5	7.5	7.5	7.5	7.5	7.5	7.5	7.5	7.5
FEL/LAT		SID 40		Table Grid		Lg Focal		AEC: Off	
cm	6	8	10	12	14	16	18	20	22
kVp	50	52	56	60	63	67	71	75	80
mAs	7.5	7.5	7.5	7.5	7.5	7.5	7.5	7.5	7.5
CAN/VD		SID 40		Table Grid		Lg Focal		AEC: Off	
cm	7	10	13	16	19	22	25	28	31
kVp	51	56	62	67	73	80	86	92	98
mAs	7.1	7.1	7.1	7.1	7.1	7.1	7.1	7.1	7.1
CAN/LAT		SID 40		Table Grid		Lg Focal		AEC: Off	
cm	7	10	13	16	19	22	25	28	31
kVp	51	56	62	67	73	80	86	92	98
mAs	7.1	7.1	7.1	7.1	7.1	7.1	7.1	7.1	7.1

Region Name: Abdomen									
FEL/LAT		SID 40		Table Grid		Lg Focal		AEC: Off	
cm	4	6	8	10	12	14	16	18	20
kVp	44	47	50	52	54	58	61	64	70
mAs	7.5	7.5	7.5	7.5	7.5	7.5	7.5	7.5	7.5
FEL/VD/DV		SID 40		Table Grid		Lg Focal		AEC: Off	
cm	4	6	8	10	12	14	16	18	20
kVp	44	47	50	52	54	58	61	64	70
mAs	7.5	7.5	7.5	7.5	7.5	7.5	7.5	7.5	7.5
CAN/LAT		SID 40		Table Grid		Lg Focal		AEC: Off	
cm	4	7	10	13	16	19	22	25	28
kVp	44	48	52	56	61	68	74	80	86
mAs	7.1	7.1	7.1	7.1	7.1	7.1	7.1	7.1	7.1
CAN/VD/DV		SID 40		Table Grid		Lg Focal		AEC: Off	
cm	4	7	10	13	16	19	22	25	28
kVp	44	48	52	56	61	68	74	80	86
mAs	7.1	7.1	7.1	7.1	7.1	7.1	7.1	7.1	7.1

