



HealthCare Equipment & Service

800-352-3201

xraysupply@yahoo.com

P.O.Box 424 Kewaskum, WI. 53040

Fax 262-251-3270

Universal Singlephase w/Digital

Radiological Position:

Cervical -AP / LAT

kV Level: 80-85

mA Level: 200

Radiological Position:

Lumbar-AP

kV Level: 100-110

mA Level: 200

Radiological Position:

Lumbar - Lat

kV Level 110-115

mA Level: 200

Width	mAs
4	2
5	2
6	3
7	3
8	5
9	6
10	7
11	7
12	10
13	10
14	10
15	12
16	15
17	18
18	18
19	20
20	30
21	30
22	35
23	40
24	40
25	50
26	60
27	75

Width	mAs
5	1
6	2
7	2
8	2
9	2
10	3
11	3
12	5
13	6
14	7
15	10
16	10
17	12
18	18
19	20
20	25
21	30
22	35
23	40
24	50
25	60
26	75
27	75
28	100
29	120

Width	mAs
4	2
5	2
6	2
7	3
8	5
9	5
10	7
11	7
12	10
13	10
14	10
15	12
16	15
17	15
18	18
19	20
20	20
21	25
22	30
23	30
24	35
25	35
26	40
27	40
28	50